














# eakin Wound Pouches™

Size\* /  
Max Volume:

## Bung Closure



45 x 30 mm — 190ml		<b>Small</b> Pack size: 10 <b>Product Code:</b> 839260	
86 x 60 mm — 180ml		<b>Small Plus</b> Pack size: 10 <b>Product Code:</b> 839267**	
110 x 75 mm — 670ml		<b>Medium</b> Pack size: 10 <b>Product Code:</b> 839261**	
175 x 110 mm — 1190ml		<b>Large</b> Pack size: 10 <b>Product Code:</b> 839262*	
245 x 160 mm — 3240ml		<b>XL Horizontal</b> Pack size: 5 <b>Product Code:</b> 839263*	<b>XL Horizontal with Access Window™</b> Pack size: 5 <b>Product Code:</b> 839264
245 x 160 mm — 3260ml		<b>XL Vertical</b> Pack size: 5 <b>Product Code:</b> 839265*	<b>XL Vertical with Access Window™</b> Pack size: 5 <b>Product Code:</b> 839268
250 x 80 mm — 1100ml		<b>Midi</b> Pack size: 5 <b>Product Code:</b> 839270*	
290 x 130 mm — 1780ml		<b>XXL</b> Pack size: 5 <b>Product Code:</b> 839266*	
250 x 385 mm — 8250ml		<b>Maxi</b> Pack size: 1 <b>Product Code:</b> 839269*	

## Fold & Tuck Closure

45 x 30 mm — 190ml		<b>Small</b> Pack size: 10 <b>Product Code:</b> 839250
110 x 75 mm — 670ml		<b>Medium</b> Pack size: 10 <b>Product Code:</b> 839251**
175 x 110 mm — 1190ml		<b>Large</b> Pack size: 10 <b>Product Code:</b> 839252*
245 x 160 mm — 3240ml		<b>XL Horizontal</b> Pack size: 5 <b>Product Code:</b> 839253*






Size:

## eakin Cohesive® skin barriers

100 x 200mm		<b>Large</b> Pack size: 5 <b>Product Code:</b> 839003
100 x 100mm		<b>Small</b> Pack size: 5 <b>Product Code:</b> 839004

Size:

## eakin® accessories

125 x 75mm		<b>eakin Access Window™ large</b> Pack size: 5 <b>Product Code:</b> 839280*	
60 x 60mm		<b>eakin Access Window™ small</b> Pack size: 5 <b>Product Code:</b> 839281**	
60g		<b>eakin Cohesive® paste</b> Pack size: 1 <b>Product Code:</b> 839010	
50ml		<b>eakin release™</b> <b>Spray, 50ml</b> Pack size: 1 <b>Product Code:</b> 839023	<b>Wipes</b> Pack size: 30 <b>Product Code:</b> 839024
60 sachets		<b>eakin perform™</b> Pack size: 60 <b>Product Code:</b> 839020	

Diameter/  
Depth:

## eakin® seals

98mm/ 2.1mm		<b>eakin Cohesive® large</b> Pack size: 10 <b>Product Code:</b> 839001
48mm/ 4.2mm		<b>eakin Cohesive® small</b> Pack size: 30 <b>Product Code:</b> 839002
48mm/ 2.8mm		<b>eakin Cohesive® slims</b> Pack size: 30 <b>Product Code:</b> 839005
85mm/ 3.0mm		<b>eakin Cohesive® stomawrap</b> Pack size: 10 <b>Product Code:</b> 839006

+ For wounds under this size.

\*Large eakin Access Window™ - can be applied to this product

\*\*Small eakin Access Window™ - can be applied to this product

15 Ballystockart Road, Comber,  
Co Down, BT23 5QY, Northern Ireland  
+44 (0)28 9187 1000 | [www.eakin.eu](http://www.eakin.eu)



© 2025 TG Eakin Ltd. | Rev 2: Jun 2024 | EWP-WPP.1 (EN)